

MS Table 2B Hi-Flow

品番: 10-1028



説明

MSテーブルは、Triple Quad™、QTRAP™、またはIVD MSシステムを収容するSciex専用製品として設計されています。MSテーブルは標準およびプレミアム形式で提供されます。プレミアムバージョンには、組み込みのInfinity 1031が含まれており、カーテン、ソース、および排気ガスの供給が可能です。さらに、MSテーブルは、これらのアプリケーションに必要な真空ポンプ用のアクセス可能な防音エンクロージャーとしても機能します。




Generator Details

- ガスの種類: 空気
- 統合型コンプレッサ: 番号

主な機能

- Sciex質量分析計を保持するために設計されたテーブル
- 単体で、または統合されたInfinity 1031コンプレッサなしの窒素ジェネレーターで供給可能
- ラボの設置面積の削減
- 最小限の設定が必要
- 防音対策済み、発電機の騒音を低減



技術仕様

MS Table 2B Hi-Flow	
ガスの種類	空気
最低入口圧力:	120psi/8.3bar
最高入口圧力:	145psi/10bar
周波数:	50 / 60 Hz
大きさ (高さx幅x奥行き) ミリメートル	804 x 1000 x 830 mm
大きさ (高さx幅x奥行き) インチ	31.7 x 39.4 x 32.7 inches

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001

